

Nur für Forschungszwecke

PHOX2A Polyklonaler Antikörper

Katalog-Nr.: 25804-1-AP



Allgemeine Informationen

Katalog-Nr.:
25804-1-AP

Größe:
150ul, Konzentration: 550 µg/ml von
Nanodrop und 300 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG22873

GenBank-Zugangsnummer:
BC041564

GeneID (NCBI):
401

Vollständiger Name:
paired-like homeobox 2a

Berechnete Masse:
284 aa, 30 kDa

Beobachtete Masse:
35-40 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB: SH-SY5Y-Zellen, HeLa-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 25804-1-AP (PHOX2A Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.