

Allgemeine Informationen

Katalog-Nr.: 25815-1-AP	GenBank-Zugangsnummer: BC031564	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 500 µg/ml von Nanodrop und 293 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 56935	Empfohlene Verdünnungen: WB 1:200-1:1000 IHC 1:50-1:500
Wirt: Kaninchen	Vollständiger Name: chromosome 11 open reading frame 75	
Isotyp: IgG	Berechnete Masse: 59 aa, 7 kDa	
Immunogen Katalognummer: AG22806	Beobachtete Masse: 7 kDa	

Anwendungen

Geprüfte Anwendungen: IHC, WB,ELISA	Positivkontrollen: WB : HEK-293-Zellen, IHC : humanes Kolonkarzinomgewebe, humanes Nierengewebe
Getestete Reaktivität: Human	
Hinweis-IHC: Antigenmaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigenmaskierung auch mit Citratpuffer pH 6,0 erfolgen.	

Hintergrundinformationen

SMCO4, also named as C11orf75 and FN5, is a Single-pass membrane and coiled-coil domain-containing protein.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




HEK-293 cells were subjected to SDS PAGE followed by western blot with 25815-1-AP (SMCO4 antibody at dilution of 1:300 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 25815-1-AP (SMCO4 antibody at dilution of 1:200 (under 40x lens).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 25815-1-AP (SMCO4 antibody at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 25815-1-AP (SMCO4 antibody at dilution of 1:200 (under 40x lens).