

Nur für Forschungszwecke

# Whirlin Polyklonaler Antikörper

Katalog-Nr.:25881-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
25881-1-AP

Größe:  
150ul, Konzentration: 800 µg/ml von  
Nanodrop und 407 µg/ml durch die  
Bradford-Methode mit BSA als  
Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG22807

GenBank-Zugangsnummer:  
BC014524

GeneID (NCBI):  
25861

Vollständiger Name:  
deafness, autosomal recessive 31

Berechnete Masse:  
907 aa, 97 kDa

Beobachtete Masse:  
46-50 kDa

Reinigungsmethode:  
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:  
WB 1:500-1:2000

## Anwendungen

Geprüfte Anwendungen:  
WB,ELISA

Getestete Reaktivität:  
Human

Positivkontrollen:

WB : COLO 320-Zellen,

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



COLO 320 cells were subjected to SDS PAGE followed by western blot with 25881-1-AP (Whirlin antibody at dilution of 1:1000 incubated at room temperature for 1.5 hours.