

Nur für Forschungszwecke

TRIM36 Polyklonaler Antikörper

Katalog-Nr.:25913-1-AP



Allgemeine Informationen

Katalog-Nr.:
25913-1-AP

Größe:

150ul, Konzentration: 1000 µg/ml von55521

Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG23240

GenBank-Zugangsnummer:

BC130334

GeneID (NCBI):

55521

Vollständiger Name:

tripartite motif-containing 36

Beobachtete Masse:

83 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : Maushodengewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse testis tissue was subjected to SDS PAGE followed by western blot with 25913-1-AP (TRIM36 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.