

Nur für Forschungszwecke

C15orf27 Polyklonaler Antikörper

Katalog-Nr.:26094-1-AP



Allgemeine Informationen

Katalog-Nr.:

26094-1-AP

Größe:

150ul, Konzentration: 450 µg/ml von Nanodrop und 320 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG23390

GenBank-Zugangsnummer:

BC104953

GeneID (NCBI):

123591

Vollständiger Name:

chromosome 15 open reading frame 27

Beobachtete Masse:

58 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB: Mausebergewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse liver tissue were subjected to SDS PAGE followed by western blot with 26094-1-AP (C15orf27 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.