

Nur für Forschungszwecke

# ECHDC2 Polyklonaler Antikörper

Katalog-Nr.:26126-1-AP



## Allgemeine Informationen

Katalog-Nr.:

26126-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG23510

GenBank-Zugangsnummer:

BC044574

GeneID (NCBI):

55268

Vollständiger Name:

enoyl Coenzyme A hydratase domain  
containing 2

Reinigungsmethode:

Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Human, Maus

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Ausgewählte Validierungsdaten



mouse liver mitochondria tissue were subjected to SDS PAGE followed by western blot with 26126-1-AP (ECHDC2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.