

Nur für Forschungszwecke

ZNF470 Polyklonaler Antikörper

Katalog-Nr.: 26134-1-AP



Allgemeine Informationen

Katalog-Nr.:

26134-1-AP

Größe:

150ul, Konzentration: 800 µg/ml von Nanodrop und 413 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG23450

GenBank-Zugangsnummer:

BC136762

GeneID (NCBI):

388566

Vollständiger Name:

zinc finger protein 470

Beobachtete Masse:

73-83 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Ratte

Positivkontrollen:

WB: Rattenhodengewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



rat testis tissue were subjected to SDS PAGE followed by western blot with 26134-1-AP (ZNF470 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.