

Nur für Forschungszwecke

RPS6KB2 Polyklonaler Antikörper

Katalog-Nr.:26194-1-AP



Allgemeine Informationen

Katalog-Nr.:
26194-1-AP

Größe:
150ul, Konzentration: 750 µg/ml von
Nanodrop und 373 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG23800

GenBank-Zugangsnummer:
BC000094

GeneID (NCBI):
6199

Vollständiger Name:
ribosomal protein S6 kinase, 70kDa,
polypeptide 2

Berechnete Masse:
53 kDa

Beobachtete Masse:
54 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human

Positivkontrollen:

WB : MCF-7-Zellen,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



MCF-7 cells were subjected to SDS PAGE followed by western blot with 26194-1-AP (RPS6KB2 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.