

Nur für Forschungszwecke

# FAM126A Polyklonaler Antikörper

Katalog-Nr.:26243-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
26243-1-AP

Größe:  
150ul, Konzentration: 1000 µg/ml von 84668  
Nanodrop und 433 µg/ml durch die  
Bradford-Methode mit BSA als  
Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG24004

GenBank-Zugangsnummer:  
BC018710

GeneID (NCBI):  
Vollständiger Name:  
family with sequence similarity 126,  
member A

Beobachtete Masse:  
50 kDa

Reinigungsmethode:  
Antigen-Affinitätsreinigung  
Empfohlene Verdünnungen:  
WB 1:500-1:2000

## Anwendungen

Geprüfte Anwendungen:  
WB, ELISA

Getestete Reaktivität:  
Human

Positivkontrollen:

WB: K-562-Zellen, A431-Zellen

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



K-562 cells were subjected to SDS PAGE followed by western blot with 26243-1-AP (FAM126A Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.