

Nur für Forschungszwecke

ANKLE2 Polyklonaler Antikörper

Katalog-Nr.:26293-1-AP



Allgemeine Informationen

Katalog-Nr.:
26293-1-AP

Größe:

150ul, Konzentration: 1000 µg/ml von23141

Nanodrop und 333 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG23387

GenBank-Zugangsnummer:

BC043157

GeneID (NCBI):

23141

Vollständiger Name:

ankyrin repeat and LEM domain
containing 2

Beobachtete Masse:

105-110 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : HeLa-Zellen, HEK-293-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



HeLa cells were subjected to SDS PAGE followed by western blot with 26293-1-AP (ANKLE2 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.