

Nur für Forschungszwecke

C10orf46 Polyklonaler Antikörper

Katalog-Nr.: 26321-1-AP



Allgemeine Informationen

Katalog-Nr.:

26321-1-AP

Größe:

150ul, Konzentration: 900 µg/ml von Nanodrop und 460 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG23694

GenBank-Zugangsnummer:

BC038294

GeneID (NCBI):

143384

Vollständiger Name:

chromosome 10 open reading frame 46

Beobachtete Masse:

41 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB: Mausebergewebe, L02-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse liver tissue were subjected to SDS PAGE followed by western blot with 26321-1-AP (C10orf46 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.