

Allgemeine Informationen

Katalog-Nr.:
26362-1-AP

Größe:
150ul, Konzentration: 650 µg/ml von
Nanodrop und 367 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG23803

GenBank-Zugangsnummer:
BC029532

GeneID (NCBI):
134187

Vollständiger Name:
POU domain class 5, transcription
factor 2

Beobachtete Masse:
42-45 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung
Empfohlene Verdünnungen:
WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human

Positivkontrollen:

WB: NCCIT cell,

Hintergrundinformationen

POU5F2, also named as POU domain, class 5, transcription factor 2m, is a 328 amino acid protein, which contains 1 homeobox DNA-binding domain and 1 POU-specific domain and belongs to the POU transcription factor family. POU5F2 as a transcription factor, localizes in the nucleus and may exert a regulatory function in meiotic events that are required for terminal differentiation of male germ cell.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



NCCIT cell were subjected to SDS PAGE followed by western blot with 26362-1-AP (POU5F2 Antibody) at dilution of 1:200 incubated at room temperature for 1.5 hours.