

Nur für Forschungszwecke

Alpha 1 microglobulin Polyklonaler Antikörper



Katalog-Nr.: 26716-1-AP

Allgemeine Informationen

Katalog-Nr.: 26716-1-AP	GenBank-Zugangsnummer: BC041593	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 1000 µg/ml von259	GeneID (NCBI): Vollständiger Name: alpha-1-microglobulin/bikunin precursor	Empfohlene Verdünnungen: WB 1:5000-1:50000
Nanodrop;	Berechnete Masse: 352 aa, 39 kDa	
Wirt: Kaninchen	Beobachtete Masse: 30 kDa	
Isotyp: IgG		
Immunogen Katalognummer: AG25239		

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : human plasma, human urine
Getestete Reaktivität: Human	

Hintergrundinformationen

Alpha 1 microglobulin, also known as HI30, is a 27-30 kDa glycoprotein, present in various body fluids. It belongs to the lipocalin superfamily of hydrophobic ligand-binding proteins forming an internal ligand-binding pocket. The protein acts as a mediator of bacterial adhesion to polymer surfaces and is involved in inhibiting renal lithogenesis.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 26716-1-AP (Alpha 1 microglobulin antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.