

Nur für Forschungszwecke

# ENTPD5 Polyklonaler Antikörper

Katalog-Nr.:26746-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
26746-1-AP

Größe:

150ul, Konzentration: 1000 µg/ml von 957

Nanodrop und 333 µg/ml durch die  
Bradford-Methode mit BSA als  
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG25022

GenBank-Zugangsnummer:

BC130485

GeneID (NCBI):

Vollständiger Name:

ectonucleoside triphosphate  
diphosphohydrolase 5

Beobachtete Masse:

47 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

IP 0.5-4.0 µg für IP und 1:200-1:1000  
für WB

IHC 1:50-1:500

## Anwendungen

Geprüfte Anwendungen:

IHC, IP, WB, ELISA

Getestete Reaktivität:

Human, Maus

**Hinweis-IHC: Antigendemaskierung mit TE-  
Puffer pH 9,0 empfohlen. (\*) Wahlweise  
kann die Antigendemaskierung auch mit  
Citratpuffer pH 6,0 erfolgen.**

Positivkontrollen:

WB: Mauslebergewebe,

IP: Mauslebergewebe,

IHC: humanes Lebergewebe,

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.


## Ausgewählte Validierungsdaten




mouse liver tissue were subjected to SDS PAGE followed by western blot with 26746-1-AP (ENTPD5 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



IP Result of anti-ENTPD5 (IP:26746-1-AP, 4ug; Detection:26746-1-AP 1:300) with mouse liver tissue lysate 6000ug.



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 26746-1-AP (ENTPD5 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 26746-1-AP (ENTPD5 Antibody) at dilution of 1:200 (under 40x lens).