

Nur für Forschungszwecke

Myogenin Polyklonaler Antikörper

Katalog-Nr.: 26762-1-AP



Allgemeine Informationen

Katalog-Nr.:

26762-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG25081

GenBank-Zugangsnummer:

BC053899

GeneID (NCBI):

4656

Vollständiger Name:

myogenin (myogenic factor 4)

Berechnete Masse:

25 kDa

Beobachtete Masse:

30 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : Maus-Skelettmuskelgewebe, Ratten-Skelettmuskelgewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 26762-1-AP (Myogenin antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.