

Nur für Forschungszwecke

APOC2 Polyklonaler Antikörper

Katalog-Nr.: 27045-1-AP



Allgemeine Informationen

Katalog-Nr.:
27045-1-AP

Größe:
150ul, Konzentration: 1000 µg/ml von 344
Nanodrop und 667 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG24650

GenBank-Zugangsnummer:
BC005348

GeneID (NCBI):

Vollständiger Name:
apolipoprotein C-II

Beobachtete Masse:
9 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human

Positivkontrollen:

WB: humanes Plasmagewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



human plasma were subjected to SDS PAGE followed by western blot with 27045-1-AP (APOC2 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.