

Nur für Forschungszwecke

C17orf49 Polyklonaler Antikörper

Katalog-Nr.: 27114-1-AP



Allgemeine Informationen

Katalog-Nr.:
27114-1-AP

Größe:

150ul, Konzentration: 1100 µg/ml von 124944

Nanodrop und 533 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG25882

GenBank-Zugangsnummer:

BC040036

GeneID (NCBI):

124944

Vollständiger Name:

chromosome 17 open reading frame
49

Beobachtete Masse:

22-25 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:3000

IP 0.5-4.0 µg für IP und 1:500-1:3000
für WB

IHC 1:50-1:500

Anwendungen

Geprüfte Anwendungen:

IHC, IP, WB, ELISA

Getestete Reaktivität:

Human, Maus

**Hinweis-IHC: Antigendemaskierung mit TE-
Puffer pH 9,0 empfohlen. (*) Wahlweise
kann die Antigendemaskierung auch mit
Citratpuffer pH 6,0 erfolgen.**

Positivkontrollen:

WB: Jurkat-Zellen, Maus-Eierstockgewebe, PC-3-
Zellen

IP: PC-3-Zellen,

IHC: humanes Kolongewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.


Ausgewählte Validierungsdaten




Jurkat cells were subjected to SDS PAGE followed by western blot with 27114-1-AP (C17orf49 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



IP result of anti-C17orf49 (IP:27114-1-AP, 4ug; Detection:27114-1-AP 1:1500) with PC-3 cells lysate 3840 ug.



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 27114-1-AP (C17orf49 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 27114-1-AP (C17orf49 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).