

Nur für Forschungszwecke

MED9 Polyklonaler Antikörper

Katalog-Nr.: 27452-1-AP



Allgemeine Informationen

Katalog-Nr.:
27452-1-AP

Größe:
150ul, Konzentration: 300 µg/ml von
Nanodrop und 233 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG26158

GenBank-Zugangsnummer:
BC000647

GeneID (NCBI):
55090

Vollständiger Name:
mediator complex subunit 9

Berechnete Masse:
146 aa, 16 kDa

Beobachtete Masse:
16 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB: HEK-293T-Zellen,

Hintergrundinformationen

MED9, also named as MED25, is a 146 amino acid protein, which belongs to the Mediator complex subunit 9 family. MED9 is a component of the Mediator complex and a coactivator that is involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. MED9 as a mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HEK-293T cells were subjected to SDS PAGE followed by western blot with 27452-1-AP (MED9 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.