

Nur für Forschungszwecke

DDX41 Polyklonaler Antikörper

Katalog-Nr.: 27500-1-AP



Allgemeine Informationen

Katalog-Nr.:
27500-1-AP

Größe:

150ul, Konzentration: 550 µg/ml von
Nanodrop und 333 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:
AG26593

GenBank-Zugangsnummer:
BC015476

GeneID (NCBI):

51428

Vollständiger Name:

DEAD (Asp-Glu-Ala-Asp) box
polypeptide 41

Berechnete Masse:

70 kDa

Beobachtete Masse:

70 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:750-1:3000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB: MCF-7-Zellen, RAW 264.7-Zellen

Hintergrundinformationen

DDX41, also named as DEAD box protein abstract homolog, is a 622 amino acid protein, which is a probable ATP-dependent RNA helicase. DDX41 is required during post-transcriptional gene expression and may be involved in pre-mRNA splicing.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



MCF-7 cells were subjected to SDS PAGE followed by western blot with 27500-1-AP (DDX41 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.