

Nur für Forschungszwecke

ANKRD17 Polyklonaler Antikörper

Katalog-Nr.:27701-1-AP



Allgemeine Informationen

Katalog-Nr.: 27701-1-AP	GenBank-Zugangsnummer: NM_001286771	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 800 µg/ml von Nanodrop;	GeneID (NCBI): 26057	Empfohlene Verdünnungen: WB 1:500-1:3000
Wirt: Kaninchen	Vollständiger Name: ankyrin repeat domain 17	
Isotyp: IgG	Berechnete Masse: 274 kDa	
Immunogen Katalognummer: AG26767	Beobachtete Masse: 270-300 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : HEK-293T-Zellen, HeLa-Zellen
Getestete Reaktivität: Human	

Hintergrundinformationen

ANKRD17 belongs to the family of ankyrin repeat-containing proteins, and contains two distinct arrays of ankyrin repeats in its amino-terminal region. It also contains a nuclear export signal, nuclear localization signal, and a cyclin-binding RXL motif. Localization of this protein to the nucleus has been shown experimentally, and interactions between this protein and cyclin-dependent kinase 2 have been observed.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 27701-1-AP (ANKRD17 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.