

Nur für Forschungszwecke

ABCB4 Polyklonaler Antikörper

Katalog-Nr.: 27726-1-AP



Allgemeine Informationen

Katalog-Nr.: 27726-1-AP	GenBank-Zugangsnummer: BC042531	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 500 µg/ml von Nanodrop;	GeneID (NCBI): 5244	Empfohlene Verdünnungen: WB 1:1000-1:5000
Wirt: Kaninchen	Vollständiger Name: ATP-binding cassette, sub-family B (MDR/TAP), member 4	
Isotyp: IgG	Berechnete Masse: 1286 aa, 142 kDa	
Immunogen Katalognummer: AG26840	Beobachtete Masse: 130-140 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : HeLa-Zellen, HepG2-Zellen, MCF-7-Zellen
Getestete Reaktivität: Human	

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 27726-1-AP (ABCB4 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.