

Nur für Forschungszwecke

PNPLA7 Polyklonaler Antikörper

Katalog-Nr.: 27787-1-AP

1 Publikationen



Allgemeine Informationen

Katalog-Nr.:

27787-1-AP

Größe:

150ul, Konzentration: 300 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG27184

GenBank-Zugangsnummer:

NM_001098537

GeneID (NCBI):

375775

Vollständiger Name:

patatin-like phospholipase domain
containing 7

Berechnete Masse:

146 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human

Zitierte Arten:

Maus

Hintergrundinformationen

PNPLA7 (Patatin-like phospholipase domain-containing protein 7), also termed NTE-related 1 (NTE-R1) or NRE, is a member of the PNPLAs family, which is conserved protein in mice, rats and humans. It serves key roles in triglyceride hydrolysis, energy metabolism, lipid droplet (LD), and in regulation of adipocyte differentiation. It is reported that PNPLA7 was down-regulated in HCC through hypermethylation of its promoter, which indicates that PNPLA7 may be a potential biomarker for HCC diagnosis (PMID: 16799181, 27347198).

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Yuki Tamura	34523264	Physiol Rep	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

