

Nur für Forschungszwecke

RASEF Polyklonaler Antikörper

Katalog-Nr.: 28162-1-AP



Allgemeine Informationen

Katalog-Nr.:
28162-1-AP

Größe:

150ul, Konzentration: 900 µg/ml von
Nanodrop und 400 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:
AG28049

GenBank-Zugangsnummer:

BC023566

GeneID (NCBI):

158158

Vollständiger Name:

RAS and EF-hand domain containing

Berechnete Masse:

740 aa, 83 kDa

Beobachtete Masse:

90-100 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:1000-1:3000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Ratte

Positivkontrollen:

WB: LNCaP-Zellen, PC-12-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 28162-1-AP (RASEF antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.