

Nur für Forschungszwecke

BPIFA2 Polyklonaler Antikörper

Katalog-Nr.:28293-1-AP



Allgemeine Informationen

Katalog-Nr.:
28293-1-AP

Größe:

150ul, Konzentration: 800 µg/ml von
Nanodrop und 367 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:
AG27401

GenBank-Zugangsnummer:

BC065726

GeneID (NCBI):

140683

Vollständiger Name:

chromosome 20 open reading frame
70

Berechnete Masse:

249 aa, 27 kDa

Beobachtete Masse:

30-37 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:5000-1:50000

Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : human saliva,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



human saliva tissue were subjected to SDS PAGE followed by western blot with 28293-1-AP (BPIFA2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.