

Nur für Forschungszwecke

CDS2 Polyklonaler Antikörper

Katalog-Nr.:28407-1-AP



Allgemeine Informationen

Katalog-Nr.:
28407-1-AP

Größe:

150ul, Konzentration: 1000 µg/ml von 8760

Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG29228

GenBank-Zugangsnummer:

BC025751

GeneID (NCBI):

8760

Vollständiger Name:

CDP-diacylglycerol synthase
(phosphatidate cytidyltransferase)

2

Berechnete Masse:

445 aa, 51 kDa

Beobachtete Masse:

50-51 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB: HeLa-Zellen, HepG2-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 28407-1-AP (CDS2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.