

Nur für Forschungszwecke

FOXJ1 Polyklonaler Antikörper

Katalog-Nr.:28443-1-AP



Allgemeine Informationen

Katalog-Nr.:
28443-1-AP

Größe:
150ul, Konzentration: 900 µg/ml von
Nanodrop;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG28294

GenBank-Zugangsnummer:
BC046460

GeneID (NCBI):
2302

Vollständiger Name:
forkhead box J1

Berechnete Masse:
45 kDa

Beobachtete Masse:
45-55 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:1000-1:5000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human

Positivkontrollen:

WB : Caco-2-Zellen, HeLa-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 28443-1-AP (FOXJ1 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.