

Nur für Forschungszwecke

# NMNAT1 Polyklonaler Antikörper

Katalog-Nr.:28493-1-AP



## Allgemeine Informationen

<b>Katalog-Nr.:</b> 28493-1-AP	<b>GenBank-Zugangsnummer:</b> BC014943	<b>Reinigungsmethode:</b> Antigen-Affinitätsreinigung
<b>Größe:</b> 150ul , Konzentration: 600 µg/ml von Nanodrop;	<b>GeneID (NCBI):</b> 64802	<b>Empfohlene Verdünnungen:</b> WB 1:500-1:2000
<b>Wirt:</b> Kaninchen	<b>Vollständiger Name:</b> nicotinamide nucleotide adenylyltransferase 1	
<b>Isotyp:</b> IgG	<b>Berechnete Masse:</b> 33 kDa	
<b>Immunogen Katalognummer:</b> AG29584	<b>Beobachtete Masse:</b> 35 kDa	

## Anwendungen

<b>Geprüfte Anwendungen:</b> WB, ELISA	<b>Positivkontrollen:</b> WB : Maus-Skelettmuskelgewebe,
<b>Getestete Reaktivität:</b> Human, Maus	

## Hintergrundinformationen

NMNAT1 is a member of the nicotinamide-nucleotide adenylyltransferases (NMNATs) which catalyze nicotinamide adenine dinucleotide (NAD) synthesis (PMID: 28445802). NMNAT is a central enzyme in NAD biosynthesis, catalyzing the formation of NAD<sup>+</sup> from nicotinamide mononucleotide (NMN) and ATP (PMID: 17402747). NMNAT1 is widely expressed with the highest levels in skeletal muscle, heart, and kidney (PMID: 11027696). Mutations in NMNAT1 have been shown associated with the LCA9 form of the retinal degeneration pathology Leber's congenital amaurosis (PMID: 22842229, 22842230).

## Lagerung

**Lagerungsbedingungen:**  
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil  
**Lagerungspuffer:**  
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.  
Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 28493-1-AP (NMNAT1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.