

RSPO4 Polyklonaler Antikörper

Katalog-Nr.: 28559-1-AP

Allgemeine Informationen

Katalog-Nr.:
28559-1-AP

Größe:

150ul, Konzentration: 750 µg/ml von
Nanodrop und 380 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG29308

GenBank-Zugangsnummer:

NM_001029871

GeneID (NCBI):

343637

Vollständiger Name:

R-spondin family, member 4

Berechnete Masse:

26 kDa

Beobachtete Masse:

~19 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:3000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB: Rattenhirngewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Rat brain tissue lysates were subjected to SDS PAGE followed by western blot with 28559-1-AP (RSP04 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.