

Nur für Forschungszwecke

# Ajuba Polyklonaler Antikörper

Katalog-Nr.:28686-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
28686-1-AP

Größe:

150ul, Konzentration: 1000 µg/ml von 84962

Nanodrop und 400 µg/ml durch die  
Bradford-Methode mit BSA als  
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG29725

GenBank-Zugangsnummer:

BC007580

GeneID (NCBI):

Vollständiger Name:

jub, ajuba homolog (Xenopus laevis)

Berechnete Masse:

57 kDa

Beobachtete Masse:

57 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:1000-1:6000

## Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB: HeLa-Zellen, HepG2-Zellen

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.

## Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 28686-1-AP (Ajuba antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.