

Nur für Forschungszwecke

SFRS6 Polyklonaler Antikörper

Katalog-Nr.: 28715-1-AP



Allgemeine Informationen

Katalog-Nr.: 28715-1-AP	GenBank-Zugangsnummer: BC006832	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 550 µg/ml von Nanodrop;	GeneID (NCBI): 6431	Empfohlene Verdünnungen: WB 1:500-1:1000
Wirt: Kaninchen	Vollständiger Name: splicing factor, arginine/serine-rich 6	
Isotyp: IgG	Berechnete Masse: 344 aa, 40 kDa	
Immunogen Katalognummer: AG28912	Beobachtete Masse: 38 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB: U2OS-Zellen,
Getestete Reaktivität: Human	

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 28715-1-AP (SFRS6 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.