

Nur für Forschungszwecke

DHX30 Polyklonaler Antikörper

Katalog-Nr.:28814-1-AP



Allgemeine Informationen

Katalog-Nr.:
28814-1-AP

Größe:
150ul, Konzentration: 500 µg/ml von
Nanodrop und 300 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG29885

GenBank-Zugangsnummer:
BC014237

GeneID (NCBI):
22907

Vollständiger Name:
DEAH (Asp-Glu-Ala-His) box
polypeptide 30

Berechnete Masse:
501 aa, 56 kDa

Beobachtete Masse:
134 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : HEK-293T-Zellen, HuH-7-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 28814-1-AP (DHX30 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.