

Nur für Forschungszwecke

DKK2 Polyklonaler Antikörper

Katalog-Nr.: 29023-1-AP



Allgemeine Informationen

Katalog-Nr.: 29023-1-AP	GenBank-Zugangsnummer: BC075078	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 500 µg/ml von Nanodrop;	GeneID (NCBI): 27123	Empfohlene Verdünnungen: WB 1:500-1:2000
Wirt: Kaninchen	Vollständiger Name: dickkopf homolog 2 (Xenopus laevis)	
Isotyp: IgG	Berechnete Masse: 259 aa, 28 kDa	
Immunogen Katalognummer: AG30596	Beobachtete Masse: 28 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : Mausherzgewebe, Rattenherzgewebe
Getestete Reaktivität: Human, Maus, Ratte	

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 29023-1-AP (DKK2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.