

Nur für Forschungszwecke

AMD1 Polyklonaler Antikörper

Katalog-Nr.:29108-1-AP



Allgemeine Informationen

Katalog-Nr.: 29108-1-AP	GenBank-Zugangsnummer: BC000171	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 650 µg/ml von Nanodrop;	GeneID (NCBI): 262	Empfohlene Verdünnungen: WB 1:1000-1:6000
Wirt: Kaninchen	Vollständiger Name: adenosylmethionine decarboxylase 1	
Isotyp: IgG	Berechnete Masse: 38 kDa, 32 kDa	
Immunogen Katalognummer: AG30751	Beobachtete Masse: 30-40 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : A549-Zellen, DU 145-Zellen
Getestete Reaktivität: Human	

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 29108-1-AP (AMD1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.