

Nur für Forschungszwecke

# HOMER2 Polyklonaler Antikörper

Katalog-Nr.:29282-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
29282-1-AP

Größe:  
150ul , Konzentration: 550 µg/ml von  
Nanodrop;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG29718

GenBank-Zugangsnummer:  
BC012109

GeneID (NCBI):  
9455

Vollständiger Name:  
homer homolog 2 (Drosophila)

Berechnete Masse:  
41 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

## Anwendungen

Geprüfte Anwendungen:  
WB,ELISA

Getestete Reaktivität:  
Human

Positivkontrollen:

WB : HEK-293-Zellen, LNCaP-Zellen, MCF-7-Zellen

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 29282-1-AP (HOMER2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.