

Nur für Forschungszwecke

WDR47 Polyklonaler Antikörper

Katalog-Nr.:29364-1-AP



Allgemeine Informationen

Katalog-Nr.: 29364-1-AP	GenBank-Zugangsnummer: BC039254	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 900 µg/ml von Nanodrop;	GeneID (NCBI): 22911	Empfohlene Verdünnungen: WB 1:2000-1:16000
Wirt: Kaninchen	Vollständiger Name: WD repeat domain 47	
Isotyp: IgG	Berechnete Masse: 919 aa, 102 kDa	
Immunogen Katalognummer: AG30853	Beobachtete Masse: 100 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : HeLa-Zellen, K-562-Zellen
Getestete Reaktivität: Human	

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern.
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 29364-1-AP (WDR47 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.