

Nur für Forschungszwecke

SH3BP4 Polyklonaler Antikörper

Katalog-Nr.:29623-1-AP



Allgemeine Informationen

Katalog-Nr.:
29623-1-AP

Größe:

150ul, Konzentration: 600 µg/ml von
Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:
AG31221

GenBank-Zugangsnummer:

BC057396

GeneID (NCBI):

23677

Vollständiger Name:

SH3-domain binding protein 4

Berechnete Masse:

963 aa, 108 kDa

Beobachtete Masse:

110 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:1000-1:8000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB : HepG2-Zellen, Mauslungengewebe, MCF-7-Zellen, Rattenlungengewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 29623-1-AP (SH3BP4 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.