

Nur für Forschungszwecke

CDT1 Polyklonaler Antikörper

Katalog-Nr.:29766-1-AP



Allgemeine Informationen

Katalog-Nr.:
29766-1-AP

Größe:

150ul, Konzentration: 450 µg/ml von
Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG31758

GenBank-Zugangsnummer:

NM_030928

GeneID (NCBI):

81620

Vollständiger Name:

chromatin licensing and DNA
replication factor 1

Berechnete Masse:

60KD

Beobachtete Masse:

65-70 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:1000-1:4000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : MCF-7-Zellen, HeLa-Zellen, HepG2-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 29766-1-AP (CDT1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.