

Nur für Forschungszwecke

BAG2 Polyklonaler Antikörper

Katalog-Nr.: 29820-1-AP



Allgemeine Informationen

Katalog-Nr.: 29820-1-AP	GenBank-Zugangsnummer: NM_004282	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 750 µg/ml von Nanodrop;	GeneID (NCBI): 9532	Empfohlene Verdünnungen: WB 1:1000-1:4000
Wirt: Kaninchen	Vollständiger Name: BCL2-associated athanogene 2	
Isotyp: IgG	Berechnete Masse: 24 kDa	
Immunogen Katalognummer: AG30873	Beobachtete Masse: 25-27 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen:
Getestete Reaktivität: Human	WB : A549-Zellen, HepG2-Zellen, Maushodengewebe, Rattenhodengewebe, U-251-Zellen, U-87 MG-Zellen

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern.
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



A549 cells were subjected to SDS PAGE followed by western blot with 29820-1-AP (BAG2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.