

Nur für Forschungszwecke

ACOT7L Polyklonaler Antikörper

Katalog-Nr.: 51016-2-AP



Allgemeine Informationen

Katalog-Nr.:

51016-2-AP

Größe:

150ul, Konzentration: 900 µg/ml von

Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

AK125299

GeneID (NCBI):

344967

Vollständiger Name:

acyl-CoA thioesterase 7 pseudogene

Berechnete Masse:

28 kDa

Beobachtete Masse:

28 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : HEK-293-Zellen, HeLa-Zellen, humanes

Hirngewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 51016-2-AP (ACOT7L antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.