

Nur für Forschungszwecke

# CRISP3 Polyklonaler Antikörper

Katalog-Nr.: 51022-1-AP



## Allgemeine Informationen

Katalog-Nr.:

51022-1-AP

Größe:

150ul, Konzentration: 133 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

BC069602

GeneID (NCBI):

10321

Vollständiger Name:

cysteine-rich secretory protein 3

Berechnete Masse:

28 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

## Hintergrundinformationen

CRISP3, also known as SGP28, is a glycoprotein that belongs to a family of cysteine-rich secretory proteins (CRISPs). CRISP3 was originally discovered in human neutrophils, but transcripts are widely distributed in exocrine glands (salivary glands, pancreas, and prostate) and also found at lower levels in epididymis, ovary, thymus, and colon. CRISP3 is a potential biomarker for prostate cancer.

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



HEK-293 cells were subjected to SDS PAGE followed by western blot with 51022-1-AP (CRISP3 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.