

Nur für Forschungszwecke

CTRP6 Polyklonaler Antikörper

Katalog-Nr.: 51030-1-AP



Allgemeine Informationen

Katalog-Nr.:
51030-1-AP

Größe:
150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

GenBank-Zugangsnummer:
BC020551

GeneID (NCBI):
114904

Vollständiger Name:
C1q and tumor necrosis factor related
protein 6

Berechnete Masse:
278 aa, 31 kDa

Beobachtete Masse:
55 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung
Empfohlene Verdünnungen:
WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA
Getestete Reaktivität:
Human

Positivkontrollen:
WB: HeLa-Zellen,

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HeLa cells were subjected to SDS PAGE followed by western blot with 51030-1-AP (CTRP6 antibody) at dilution of 1:0 incubated at room temperature for 1.5 hours.