

Nur für Forschungszwecke

Cadm1 Polyklonaler Antikörper

Katalog-Nr.:51059-1-AP



Allgemeine Informationen

Katalog-Nr.:
51059-1-AP

Größe:
150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG0449

GenBank-Zugangsnummer:
AF539424

GeneID (NCBI):
54725

Vollständiger Name:
cell adhesion molecule 1

Berechnete Masse:
49 kDa

Beobachtete Masse:
49 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Maus

Positivkontrollen:

WB: Maushodengewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern.

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse testis tissue were subjected to SDS PAGE followed by western blot with 51059-1-AP (Cadm1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.