

Allgemeine Informationen

Katalog-Nr.:

51101-1-AP

Größe:

150ul, Konzentration: 500 µg/ml von Nanodrop und 333 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

BC104994

GeneID (NCBI):

23325

Vollständiger Name:

KIAA1033

Berechnete Masse:

136.4 kDa

Beobachtete Masse:

136 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:3000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human

Zitierte Arten:

Human, Maus

Positivkontrollen:

WB: A549-Zellen, HeLa-Zellen

Hintergrundinformationen

KIAA1033, also named as WASH-interacting protein and SWIP, is a component of WASH complex. The WASH complex plays a key role in the fission of tubules that serve as transport intermediates during endosome sorting. KIAA1033 has two isoforms with MW 100 kDa and 136 kDa.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Benjamin J Tyrrell	27390154	Pigment Cell Melanoma Res	WB
Yifei Pei	35993307	Autophagy	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 51101-1-AP (KIAA1033 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.