

Nur für Forschungszwecke

# CD11d/Integrin Alpha D Polyklonaler Antikörper



Katalog-Nr.:55162-1-AP

## Allgemeine Informationen

Katalog-Nr.:  
55162-1-AP

Größe:  
150ul , Konzentration: 900 µg/ml von  
Nanodrop und 467 µg/ml durch die  
Bradford-Methode mit BSA als  
Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

GenBank-Zugangsnummer:  
NM\_005353

GeneID (NCBI):  
3681

Vollständiger Name:  
integrin, alpha D

Berechnete Masse:  
127 kDa

Beobachtete Masse:  
126-150 kDa

Reinigungsmethode:  
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:  
WB 1:500-1:1000

## Anwendungen

Geprüfte Anwendungen:  
WB,ELISA

Getestete Reaktivität:  
Human

Positivkontrollen:

WB : Jurkat-Zellen,

## Hintergrundinformationen

CD11d, also named as ITGAD and ADB2, belongs to the integrin alpha chain family. It is a receptor for ICAM3 and VCAM1. CD11d may play a role in the atherosclerotic process such as clearing lipoproteins from plaques and in phagocytosis of blood-borne pathogens, particulate matter, and senescent erythrocytes from the blood. This antibody is specific to CD11d.

## Lagerung

Lagerungsbedingungen:  
Bei -20°C lagern.

Lagerungspuffer:  
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



Jurkat cells were subjected to SDS PAGE followed by western blot with 55162-1-AP (CD11d antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.