

Nur für Forschungszwecke

BNC1 Polyklonaler Antikörper

Katalog-Nr.: 55233-1-AP



Allgemeine Informationen

Katalog-Nr.:
55233-1-AP

Größe:
150ul, Konzentration: 413 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

GenBank-Zugangsnummer:
NM_001717

GeneID (NCBI):
646

Vollständiger Name:
basonuclein 1

Berechnete Masse:
111 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung
Empfohlene Verdünnungen:
IF 1:50-1:500

Anwendungen

Geprüfte Anwendungen:
IF, ELISA

Getestete Reaktivität:
Human, Maus

Positivkontrollen:
IF: HepG2-Zellen,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern.

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using 55233-1-AP (BNC1 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).