

NUMB Polyklonaler Antikörper

Katalog-Nr.: 55321-1-AP

Allgemeine Informationen

Katalog-Nr.:

55321-1-AP

Größe:

150ul, Konzentration: 233 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_003744

GeneID (NCBI):

8650

Vollständiger Name:

numb homolog (Drosophila)

Berechnete Masse:

71 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

NUMB plays a role in the process of neurogenesis. It is required throughout embryonic neurogenesis to maintain neural progenitor cells, also called radial glial cells (RGCs), by allowing their daughter cells to choose progenitor over neuronal cell fate. NUMB may also mediate local repair of brain ventricular wall damage. NUMB may play a role in further understanding the NUMB function in development and cancer. This antibody can recognize isoform 1, 2, 5 and 7.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

