

Nur für Forschungszwecke

p140Cap Polyklonaler Antikörper

Katalog-Nr.:55362-1-AP

1 Publikationen



Allgemeine Informationen

Katalog-Nr.:
55362-1-AP

Größe:

150ul , Konzentration: 1000 µg/ml von80725

Nanodrop und 433 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_025248

GeneID (NCBI):

80725

Vollständiger Name:

SNAP25-interacting protein

Berechnete Masse:

127 kDa

Beobachtete Masse:

140-145 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:300-1:1000

IP 0.5-4.0 ug für IP und 1:500-1:1000

für WB

IHC 1:50-1:500

Anwendungen

Geprüfte Anwendungen:

IHC, IP, WB, ELISA

In Publikationen genannte Anwendungen:

IHC, WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Human

**Hinweis-IHC: Antigendemaskierung mit TE-
Puffer pH 9,0 empfohlen. (*) Wahlweise
kann die Antigendemaskierung auch mit
Citratpuffer pH 6,0 erfolgen.**

Positivkontrollen:

WB : Maushirngewebe, Rattenhirngewebe

IP : Maushirngewebe,

IHC : Maushirngewebe,

Hintergrundinformationen

Bemerkenswerte Veröffentlichungen

Verfasser

Pubmed ID

Journal

Anwendung

Wen Wen Xu

37149929

EBioMedicine

WB,IHC

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

**This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.**


Ausgewählte Validierungsdaten




Various lysates were subjected to SDS PAGE followed by western blot with 55362-1-AP (p140Cap antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



IP result of anti-p140CaP (IP:55362-1-AP, 4 μ g; Detection:55362-1-AP 1:500) with mouse brain tissue lysate 4000 μ g.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 55362-1-AP (p140Cap antibody) at dilution of 1:200 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 55362-1-AP (p140Cap antibody) at dilution of 1:200 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).