

Nur für Forschungszwecke

Cdc20 Monoklonaler Antikörper

Katalog-Nr.: 60010-1-Ig



Allgemeine Informationen

Katalog-Nr.:

60010-1-Ig

Größe:

150ul, Konzentration: 1000 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Maus

Isotyp:

IgG2c

Immunogen Katalognummer:

AG0295

GenBank-Zugangsnummer:

BC001088

GeneID (NCBI):

991

Vollständiger Name:

cell division cycle 20 homolog (S.
cerevisiae)

Berechnete Masse:

55 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-
Präzipitation

CloneNo.:

1E8A8

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

