

Nur für Forschungszwecke

GATA1 Monoklonaler Antikörper

Katalog-Nr.:60011-2-Ig



Allgemeine Informationen

Katalog-Nr.:

60011-2-Ig

Größe:

150ul

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG1350

GenBank-Zugangsnummer:

BC009797

GeneID (NCBI):

2623

Vollständiger Name:

GATA binding protein 1 (globin transcription factor 1)

Berechnete Masse:

43 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-Präzipitation

CloneNo.:

5B1A7

Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

