

Nur für Forschungszwecke

IFNAR2 Monoklonaler Antikörper

Katalog-Nr.:60013-2-Ig



Allgemeine Informationen

Katalog-Nr.:

60013-2-Ig

Größe:

150ul

Wirt:

Maus

Isotyp:

IgM

Immunogen Katalognummer:

AG0191

GenBank-Zugangsnummer:

BC002793

GeneID (NCBI):

3455

Vollständiger Name:

interferon (alpha, beta and omega)

receptor 2

Berechnete Masse:

331 aa, 37 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-
Präzipitation

CloneNo.:

4G8G1

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

